

## EXHIBIT A

---

Project No. 612.2/1295Xf61

Date Started

Object: Preparation of enzyme degraded Guar and Locust Bean Gum.

Guar Gum

2 lb of Guar Gum

Solids = 20%

Sodium Benzoate = 0.1% solution

ph = 4, T = 70°C

Graminase enzyme

Locust Bean Gum

2 lb of Locust Bean Gum

Solids = 20%

Sodium Benzoate = 0.1% solution

ph = 4, T = 70°C

Graminase enzyme

Procedure: 0.1% Sodium Benzoate dissolved in deion. H<sub>2</sub>O, ph = 4, tied to flask and placed in constant temp. bath of bath. 70°C then added 10 ml of Graminase 1 lb of Guar Gum added 2 lb of Guar Gum added 10 ml Graminase 1 lb Locust Bean Gum and then added 2 lb of Locust Bean Gum and recd. H<sub>2</sub>O-25% L.S.

SAMPLE	Amount of Enzyme	Rx Time (hrs)	% PS	SAMPLE	Amount of Enzyme	Rx Time (hrs)	% PS
6231/134-1 (Guar Gum)	10 ml / 1 lb Gum	0/14	12% PS	6231/134-2 (Locust Bean Gum)	10 ml / 1 lb Gum	1 hr.	25% PS
	added 7 ml 1 lb of gum	0/45	16% PS				Used 20 ml of 1% soln. + deion. water 5:1 ratio solution - ph = 2
	added 5 ml / 1 lb gum	0/14	17% PS				used 3 ml of 1% soln. + deion. water 5:1 ratio solution - ph = 2

OBSERVATION:

These larger quantities were prepared to examine the functional properties of these materials in food systems.

WORK OF:

Jane Smith

WITNESS THIS DOCUMENT AND UNDERSTAND ITS CONTENTS

DATE:

M. Smith

Signature

Date

After witnessing, corrections or changes may never be made in the grid area, but may be noted at any time in this margin.

Project No. 631.2/1412+EE

Date Started

Object

Preparation of enzyme degraded Amarant.  
(big batch)

Materials:

2 1/2 lb of Amarant gum  
0.1% Sodium BenzateSolids  $\rightarrow$  20%  
 $\rightarrow$  58°C  
Cellulast enzyme (about 25 ml / 1 lb gum)  
ph  $\rightarrow$  5 (6231:137-1)300 g of Amarant gum (new batch)  
0.1% Sodium BenzateSolids  $\rightarrow$  20%  
 $\rightarrow$  58°C  
Cellulast enzyme (about 20 ml / 1 lb of gum)

Procedure:

Sample	Amount of enzyme used	Rx Time (hrs)	% RS
6231:137			
6231:137-1	25 ml of Cellulast enzyme / 1 lb of gum	20 hrs	$\rightarrow$ 22% degraded by HCl (3%)
6231:137-2	20 ml of Cellulast enzyme / 1 lb of gum	20 hrs	$\rightarrow$ 23% degraded by HCl (3%)

Observation: The big batch of Amarant 6231:137-1 will be submitted to Food Lab (along with gum and low-batch gum) to examine the functional properties of these products (also evaluated in cake and puddings).

Sana/ank

M. Hurley